

10-year status report Filtered Fume Hoods

David Withee Marketing Rep, GreenFumeHood Technology



Learning Objectives

- What a "filtered" fume hood is and is not
- Benefits of using filtered hoods
- Applications for filtered hoods





That's a lot of conditioned air!

- Estimated 750,000 ducted fume hoods operating in the USA
- In general, operating 24/7
- Utility cost approximately\$4.2 billion/year





So, let's slay the lab's biggest energy hog!






18 years ago . . .

- BIG change to fume hoods
- Still controversial
- Low Flow

Low Flow (High Performance) Fume Hoods

Utilize fume hoods that are designed to operate safely at lower face velocities, i.e., 70 FPM rather than 100 FPM.





10 years ago . . .

- Filtered fume hoods
- Emphasizing safety and sustainability.





DRIVER **#1**

- Pollution is pollution, regardless of dilution!
- Where does fume hood exhaust go?
- What is your future liability?




Driver # 2:

Energy savings






Driver #3: Mechanical systems





Cost per Sq.Ft. Section w/DUCTED Hood Images & Data credit: Brewer Construction, Houston, TX





Cost per Sq.Ft. Section w/FILTERED Hood Images & Data credit: Brewer Construction, Houston, TX

NOT a ductless hood

 Single filter cannot handle acids AND bases AND solvents





A *filtered* fume hood

 <u>Multiple</u> safety factors not found in standard or ductless fume hoods





ANSI Z9.5

Section 5.3.6.2 Exhaust Hood Air

➢ Reads

- ✓ Exhaust hood air shall not be re-circulated to other areas
- ✓ Hood exhaust air meeting the same criteria as noted in Section 5.3.6.1 <u>shall only be re-</u> <u>circulated to the same work area</u>




Multiple modules per hood for layers of safety





Unparalleled safety features


- Independent fan/filtration columns
- Redundant filters on top
- Air sensors between filters
- All electronics monitored multiple times per second
- If anything is amiss, alarm sounds





Detection at <1% Daily TLV Exposure

- Suite of Detectors:
 - Acid
 - Solvent
 - Lab Ambient Air
 - Temperature
 - Sash Sensor







(side view of hood)

Activation Process

- Organized structure
- Steam & Heat create spaces between carbon layers
- 15,000 ft² surface area per gram







Filters for science

- Neutrodine filtration technology retains:
 - Polar organic solvents
 - Non-polar organic solvents
 - Inorganic bases
 - Inorganic acids







Addressing Lab Sustainability



- Reduce carbon footprint
- Reduce operational costs
- Reduce construction cost and time
- Flexibility for changes and growth

While IMPROVING safety!

North America:

(550)+ GreenFumeHoods (200)+ Projects 350,000+ CFM Saved/Yr \$2+ mil Saved/Yr

Worldwide:

(1,000)+ GreenFumeHoods (550)+ Projects





--1 er it

University of Rochester, NY (13) GreenFumeHoods

University of Rochester, NY

Type: Renovation, \$1.5mil, 6 months Date: Mar-August, 2011 Application: Organic Chemistry Teaching Lab Architect: Bergmann Associates (13) GreenFumeHoods (2) Ducted Fume Hoods (dispensing) NYSERDA Rebate: \$36,600 \$136,100 Annual Operating Cost Savings 300% Increase in fume hoods



1 18



Butler University, IN

Type: Renovation Date: August, 2011 Application: (4) Organic and (2) General **Chemistry Teaching Labs** Architect: BSA LifeStructures MEP Engineer: Greenman-Pedersen, Inc. (27) GreenFumeHoods **10% Reduction in Total First Costs 3 Month Reduction in Construction Schedule** \$33,000 Annual Operating Cost Savings



Marywood University, PA (8) GreenFumeHoods

image by Ellenzweig Architects

Marywood University, PA

Type: Gut Renovation Date: October, 2012 Application: Organic Chemistry Teaching Lab Architect: Ellenzweig MEP Engineer: Greenman-Pedersen, Inc. (8) GreenFumeHoods 30% Reduction in HVAC Systems Lower Total First Costs \$12,000 Annual Operating Cost Savings

McMaster University, ON (8) GreenFumeHoods

.....

McMaster University, ON

Type: Renovation Date: August, 2013 Application: Integrated Science Teaching Lab (8) GreenFumeHoods (1) Ducted Fume Hood 70% Reduction in Operating Costs Wharton QS Stars 2014 award-winning program





Framingham State University, MA (49) GreenFumeHoods

1T

-

Framingham State University, MA

Type: Addition Date: May, 2014 Application: (8) Chemistry Teaching Labs Architect: Ellenzweig MEP Engineer: BR+A Engineers (49) GreenFumeHoods 25% Reduction in HVAC Systems \$330,000 Lower Total First Costs \$40,000 Annual Operating Cost Savings



JOHN J. SBRECA


T

Bristol Community College – ZNE Lab Building (13) GreenFumeHoods

Image credit: Edward Caruso Photography

THE PATH TO ZERO NET ENERGY

Bristol Community College, John J. Sbrega Science & Health Science Building





Bridgestone Tech. Center, OH (11) GreenFumeHoods

S O L HARRIS/DAY architecture

Bridgestone Tech. Center, OH

Type: New Construction Date: October, 2012 Application: R&D Lab Architect: SolHarris/Day Architecture MEP Engineer: Lee Good and Associates (11) GreenFumeHoods (1) Floor-mounted ducted fume hood 32% Reduction in Energy Consumption \$60,000 Annual Operating Cost Savings







Carleton College – Integrated Sciences Facility (11) 4 and 6-foot Protector Echos by Labconco

image by EYP

Hulings Second Floor / Olin First

Parketter 3.4

Type: "Swing Space" Renovation Date: Aug, 2017 (11) GreenFumeHoods by Labconco Convert Biology Labs to Chemistry Relocate Hoods to New Building in 2019

Use Legend

Building Suppor

Classesore

Detroited Support

Built Drange to Rem

Exist to Renain

Faculty/Shidert S

Office

Carleton College – Integrated Sciences Facility (11) 4 and 6-foot Protector Echos by Labconco

image by EYP

Carleton College – Integrated Sciences Facility (11) 4 and 6-foot Protector Echos by Labconco

Ech

image by Haldeman Homme

